

# Anti-Beclin 1 (4F4) Mouse antibody



Catalog Number: 168125

## Product name

Monoclonal

Anti-Beclin 1 (4F4) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse, Rat, Monkey, Hamster

## Applications

WB

## Antibody description

Beclin 1 (4F4) Mouse monoclonal antibody

## Dilutions

WB: 1/1000


## Preparation

Antigen: Purified recombinant human Beclin 1 protein fragments expressed in E.coli.

## Validations

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Beclin 1 in HeLa, 3T3, C6, COS7, CHO-K1, Rat brain and Jurkat cell lysates using Beclin 1 mouse mAb (1:1000 diluted). Predicted band size: 52KDa. Observed band size: 55KDa.

## Clonality