

SEPT8 antibody

Catalog Number: 107525

Product name

SEPT8 antibody

Specificity

Human; other species not tested.

Antibody description

SEPT8 Mouse Monoclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 50-58 kDa

Preparation

This antibody was obtained by immunization of SEPT8 recombinant protein (Accession Number: BC001329). Purification method: Thiophilic affinity chromatograph.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107525(SEPT8 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107525(SEPT8 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)