

Anti-JNK3 (3C9) Mouse antibody



Catalog Number: 168525

Product name

Anti-JNK3 (3C9) Mouse antibody

Dilutions

WB: 1/500 - 1/2000

Specificity

ICC: 1/200 - 1/1000

Human, Mouse

ELISA: 1/10000

Antibody description

JNK3 (3C9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human MAPK10 (aa28-233) expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB, ICC/IF

Validations



Immunofluorescence staining of methanol-fixed A431 (left) and Hela (right) cells showing cytoplasmic and membrane localization.



Western blot detection of JNK3 in Hela, THP-1 and A549 cell lysates using JNK3 mouse mAb (1:1000 diluted). Predicted band size: 54 kDa. Observed band size: 54 kDa.