

Anti-JNK3 (3C9) Mouse antibody



Catalog Number: 168525

Product name

Anti-JNK3 (3C9) Mouse antibody

Specificity

Human, Mouse

Antibody description

JNK3 (3C9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human MAPK10 (aa28-233) expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, ICC/IF


Dilutions

WB: 1/500 - 1/2000


ICC: 1/200 - 1/1000

ELISA: 1/10000

Validations



Immunofluorescence staining of methanol-fixed A431 (left) and HeLa (right) cells showing cytoplasmic and membrane localization.



Western blot detection of JNK3 in HeLa, THP-1 and A431 cell lysates using JNK3 mouse mAb (1:1000 diluted). Predicted band size: 54KDa. Observed band size: 54KDa.