

MYEOV2 antibody



Catalog Number: 112925

Product name

Polyclonal

MYEOV2 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

MYEOV2 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells. Observed molecular weight by Western-blot: 27kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of MYEOV2 recombinant protein (Accession Number: BC033955). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112925(MYEOV2 Antibody) at dilution of 1:300