

# Anti-N Cadherin Rabbit antibody



Catalog Number: 166925

## Product name

Rabbit IgG

Anti-N Cadherin Rabbit antibody

## Applications

## Specificity

WB, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

N Cadherin Rabbit polyclonal antibody

FC: 1/100

## Preparation

IP: 1/20

Antigen: Recombinant protein of human N Cadherin

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type



Western blot detection of N Cadherin in HepG2 cell lysates using N Cadherin Rabbit pAb(1:1000 diluted).Predicted band size:100KDa.Observed band size:140KDa.