

AK3L1 antibody

Catalog Number: 107925

Product name

AK3L1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AK3L1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, human liver tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of AK3L1 recombinant protein (Accession Number: NM_203464). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107925(AK3L1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107925(AK3L1 antibody) at dilution of 1:100 (under 10x lens)