

# AK3L1 antibody



Catalog Number: 107925

## Product name

AK3L1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AK3L1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, human liver tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of AK3L1 recombinant protein (Accession Number: NM\_203464). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107925(AK3L1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107925(AK3L1 antibody) at dilution of 1:100 (under 10x lens)