

Anti-hnRNP Q Rabbit antibody



Catalog Number: 168925

Product name

Anti-hnRNP Q Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

hnRNP Q Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human hnRNP Q

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

1:2000-1:10000


1:20-1:50

1:50

Validations



Immunocytochemistry of hnRNP Q(green) in HeLa cells using hnRNP Q Rabbit mAb at dilution 1/200, and DAPI(blue)



and DAPI(blue)

Western blot detection of hnRNP Q in K562 cell lysates using hnRNP Q Rabbit mAb(1:1000 diluted). Predicted band size: 70kDa. Observed band size: 70kDa.