

CAMKV antibody



Catalog Number: 108925

Product name

CAMKV antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CAMKV Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells. Observed molecular weight by Western-blot: 54-75kd

Preparation

This antibody was obtained by immunization of CAMKV recombinant protein (Accession Number: NM_024046). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108925(CAMKV antibody) at dilution of 1:1200