

Human ADPRH/ARH1 (His tag) recombinant protein



Catalog Number: 519925

General Information

Gene Name Synonym

[Protein ADP-Ribosylarginine] Hydrolase, ADP-Ribosylarginine Hydrolase, ADP-Ribose-L-Arginine Cleaving Enzyme, ADPRH, ARH1

Protein Construction

A DNA sequence encoding the Human ADPRH/ARH1 Met1-Leu357 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 50% Glycerol, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ADPRH/ARH1 (His tag) recombinant protein