

CACNB3 Polyclonal Antibody



Catalog Number: 162225

Product name

CACNB3 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to CACNB3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 355-484 of human CACNB3 (NP_000716.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - CACNB3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CACNB3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - CACNB3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat ovary using CACNB3 antibody at dilution of 1:150 (40x lens).