

# Anti-Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit antibody



Catalog Number: 167225

## Product name

Anti-Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Heterogeneous Nuclear Ribonucleoprotein (A1-like)

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/100

FC: 1/20

## Validations



Western blot detection of Heterogeneous Nuclear Ribonucleoprotein (A1-like) in K562, C6, 3T3, HeLa cell lysates using Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit pAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.