

Anti-Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit antibody



Catalog Number: 167225

Product name

Anti-Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Heterogeneous Nuclear Ribonucleoprotein (A1-like)

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/100

FC: 1/20

Validations



Western blot detection of Heterogeneous Nuclear Ribonucleoprotein (A1-like) in K562, C6, 3T3, HeLa cell lysates using Heterogeneous Nuclear Ribonucleoprotein (A1-like) Rabbit pAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.