

Anti-Cystatin C (6F12) Mouse antibody



Catalog Number: 168225

Product name

Anti-Cystatin C (6F12) Mouse antibody

Specificity

Human

Antibody description

Cystatin C (6F12) Mouse monoclonal antibody

Preparation

Antigen: Recombinant human cystatin c protein.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1

Applications

sELISA Capture

Validations



Standard Curve for Human Cystatin C: Capture Antibody Mouse mAb [6F12-C7-D8] to Human Cystatin C at 4u03bcg/ml and Detector Antibody Mouse mAb [7F6-A5-F3]to Human Cystatin C at 0.5u03bcg/ml.