

# ANGPTL4 Polyclonal Antibody



Catalog Number: 161625

---

## Product name

ANGPTL4 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ANGPTL4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 26-406 of human ANGPTL4 (NP\_647475.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - ANGPTL4 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ANGPTL4 antibody . Blue: DAPI for nuclear staining.