

CTBP2 antibody

Catalog Number: 109625

Product name

CTBP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CTBP2 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of CTBP2 recombinant protein (Accession Number: XM_017015756). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109625(CTBP2 antibody) at dilution of 1:800



IP Result of anti-CTBP2 (IP:Catalog No:109625, 4ug; Detection:Catalog No:109625 1:500) with mouse heart tissue lysate 3200ug.