

PLA2G5 Polyclonal Antibody



Catalog Number: 163325

Product name

PLA2G5 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to PLA2G5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human PLA2G5 (NP_000920.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - PLA2G5 Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using PLA2G5 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.