

Human RAET1E / ULBP4 (His Tag) recombinant protein



Catalog Number: 503325

General Information

Gene Name Synonym

Lymphocyte effector toxicity activation ligand; RAE-1-like transcript 4; Retinoic acid early transcript 1E

Protein Construction

A DNA sequence encoding the human RAET1E (NP_631904.1) (Met1-Asp225) was expressed with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 80 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His 31

Molecular Mass

The recombinant human RAET1E consists 206 amino acids and predicts a molecular mass of 23.8 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human RAET1E / ULBP4 Protein (His Tag) SDS-PAGE