

# Anti-CDK1 (4E2) Mouse antibody



Catalog Number: 168325

## Product name

Mouse IgG1

Anti-CDK1 (4E2) Mouse antibody

## Applications

## Specificity

WB

Human, Mouse, Rat, Monkey

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

CDK1 (4E2) Mouse monoclonal antibody

ELISA: 1/10000

## Preparation

Antigen: Purified recombinant fragment of human CDK1 expressed in E. Coli.

## Validations

## Formulation

Ascitic fluid containing 0.03% sodium azide.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of CDK1 in A431, Ramos, COS7 and HeLa cell lysates using CDK1 mouse mAb (1:500 diluted). Predicted band size: 34KDa. Observed band size: 34KDa.