

MRP4 / ABCC4 Polyclonal Antibody



Catalog Number: 163725

Product name

MRP4 / ABCC4 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MRP4 / ABCC4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 495-705 of human MRP4 / ABCC4 (NP_001098985.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - MRP4 / ABCC4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MRP4 / ABCC4 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunofluorescence - MRP4 / ABCC4 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MRP4 / ABCC4 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.