

Anti-STAT1 Rabbit antibody



Catalog Number: 166725

Product name

Anti-STAT1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

STAT1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human STAT1 alpha

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/1000

FC: 1/20

IP: 1/20

Validations



Western blot detection of STAT1 alpha in K562, Rat Brain, HeLa cell lysates using STAT1 alpha Rabbit pAb(1:1000 diluted). Predicted band size: 87kDa. Observed band size: 87kDa.