

Anti-XAF1 antibody



Catalog Number: 176725

Product name

Anti-XAF1 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to XAF1

Preparation

This antigen of this antibody was amino acids 1-225 mapping at the n-terminus of xaf1 of human origin.

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 40%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IP, IF

Dilutions

WB: 1:100-1:1,000

IP: 1-2 µg per 100-500 µg of total protein(1 ml of cell lysate)

IF: 1:50-1:500

Validations




Fig1: Western blot analysis of XAF1 expression in non-transfected (A) and human XAF1 transfected (B) 293T whole cell lysates.