

# Anti-DGKZ/DGK-zeta Rabbit antibody



Catalog Number: 167725

## Product name

Anti-DGKZ/DGK-zeta Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

DGKZ/DGK-zeta Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human DGKZ/DGK-zeta

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000

IHC: 1/500

ICC/IF: 1/20

FC: 1/50

## Validations



Western blot detection of DGKZ/DGK-zeta in HeLa, A549, HL-60, U2OS cell lysates using DGKZ/DGK-zeta Rabbit mAb(1:1000 diluted). Predicted band size: 104kDa. Observed band size: 124, 114kDa.