

ACAP2 antibody

Catalog Number: 107725

Product name

ACAP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACAP2 Rabbit Polyclonal antibody. Positive IP detected in NIH/3T3 cells. Positive WB detected in NIH/3T3 cells, Jurkat cells. Observed molecular weight by Western-blot: 88 kDa

Preparation

This antibody was obtained by immunization of ACAP2 recombinant protein (Accession Number: NM_012287). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107725(ACAP2 antibody) at dilution of 1:300



IP Result of anti-ACAP2 (IP:Catalog No:107725, 4ug; Detection:Catalog No:107725 1:300) with NIH/3T3 cells lysate 3600ug.