

# TH1L antibody

Catalog Number: 116035

---

**Product name**

TH1L antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

TH1L Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 66 kDa

**Preparation**

This antibody was obtained by immunization of TH1L recombinant protein (Accession Number: NM\_016397). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:500-1:5000

**Validations**

A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116035(TH1L antibody) at dilution of 1:500