

Anti-FUBP1/FPB Rabbit antibody



Catalog Number: 167035

Product name

Anti-FUBP1/FPB Rabbit antibody

WB: 1/500-1/2000

IHC: 1/50-1/200

Specificity

Human, Mouse, Rat

ICC/IF: 1/50-1/200

FC: 1/50

Antibody description

FUBP1/FPB Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human FUBP1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC

Dilutions

Validations



Western blot analysis of extracts from Jurkat, K562 and HeLa cells using FUBP1 Rabbit pAb at 1:1000 dilution. Predicted band size: 74kDa. Observed band size: 74kDa.



Immunohistochemical analysis of paraffin-embedded human spleen, using FUBP1 Antibody.