

# ZDHHC7 antibody



Catalog Number: 117035

## Product name

ZDHHC7 antibody

## Immunogen

[Human ZDHHC7 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ZDHHC7 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, MCF-7 cells. Observed molecular weight by Western-blot: 45-49 kDa

## Preparation

This antibody was obtained by immunization of ZDHHC7 recombinant protein (Accession Number: NM\_001145548). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:117035(ZDHHC7 Antibody) at dilution of 1:1000