

APOLD1 antibody

Catalog Number: 108035

Product name

APOLD1 antibody

Immunogen

[Human APOLD1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

APOLD1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of APOLD1 recombinant protein (Accession Number: NM_030817). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108035(APOLD1 Antibody) at dilution of 1:1000