

SMARCD3 antibody



Catalog Number: 115435

Product name

SMARCD3 antibody

Immunogen

[Human SMARCD3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SMARCD3 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of SMARCD3 recombinant protein (Accession Number: NM_003078). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115435(SMARCD3 antibody) at dilution of 1:400



Immunofluorescent analysis of HeLa cells, using SMARCD3 antibody Catalog No:115435 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG (green). Blue pseudocolor = DAPI (fluorescent DNA dye).