

TSLP antibody

Catalog Number: 116435

Product name

TSLP antibody

Specificity

human; other species not tested.

Antibody description

TSLP Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer.

Preparation

This antibody was obtained by immunization of TSLP recombinant protein (Accession Number: BC016720). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116435(TSLP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116435(TSLP antibody) at dilution of 1:50 (under 40x lens)