

# Biotin OPTN antibody



Catalog Number: 107435

## Product name

Biotin OPTN antibody

## Specificity

Human; other species not tested.

## Antibody description

Biotin OPTN Mouse Monoclonal antibody. Positive WB detected in human brain. Observed molecular weight by Western-blot: 70-74 kDa

## Preparation

This antibody was obtained by immunization of OPTN recombinant protein (Accession Number: NM\_001008212). Purification method: Protein G purified.

## Formulation

PBS(10mM), with 50% glycerol, without sodium azide.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Fetal human brain was subjected to SDS-PAGE followed by western blot with Biotin-60293 (biotin labeled OPTN monoclonal antibody) at dilution of 1:1000. Avidin-HRP was used instead of secondary antibody.