

Human DUSP3/VHR (His tag) recombinant protein



Catalog Number: 518435

General Information

Gene Name Synonym

Dual specificity protein phosphatase
3,DUSP3,Dual specificity protein phosphatase
VHR,Vaccinia H1-related phosphatase,VHR

Protein Construction

A DNA sequence encoding the Human
DUSP3/VHR Ser2-Pro185 is expressed with a 6His
tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of
PBS,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human DUSP3/VHR (His tag) recombinant protein