

Human DUSP3/VHR (His tag) recombinant protein



Catalog Number: 518435

General Information

Gene Name Synonym

Dual specificity protein phosphatase 3, DUSP3, Dual specificity protein phosphatase VHR, Vaccinia H1-related phosphatase, VHR

Protein Construction

A DNA sequence encoding the Human DUSP3/VHR Ser2-Pro185 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human DUSP3/VHR (His tag) recombinant protein