

Anti-ABCD4 antibody

Catalog Number: 175835

Product name

Anti-ABCD4 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ABCD4

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human abcd4. aa range:111-160

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using ABCD4 Polyclonal Antibody diluted at 1: 500