

Anti-MCAK Rabbit antibody



Catalog Number: 167835

Product name

Anti-MCAK Rabbit antibody

Specificity

Human

Antibody description

MCAK Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MCAK

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/20

FC: 1/50

Validations



Western blot detection of MCAK in HeLa,CHO-K1,Ramos cell lysates using MCAK Rabbit mAb(1:1000 diluted).Predicted band size:81kDa.Observed band size:81kDa.