

# Human PPP1R1A/IPP1 (His tag) recombinant protein



Catalog Number: 518835

## General Information

### Gene Name Synonym

Protein Phosphatase 1 Regulatory Subunit 1A, Protein Phosphatase Inhibitor 1, I-1, IPP-1, PPP1R1A, IPP1

### Protein Construction

A DNA sequence encoding the Human PPP1R1A/IPP1 Met1-Val171 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 50% Glycerol, pH 8.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PPP1R1A/IPP1 (His tag) recombinant protein