

DDX5,p68 antibody



Catalog Number: 109835

Product name

DDX5,p68 antibody

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

DDX5,p68 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in NIH/3T3 cells, A431 cells, HeLa cells. Observed molecular weight by Western-blot: 65-68kd

Preparation

This antibody was obtained by immunization of DDX5,p68 recombinant protein (Accession Number: NM_001320596). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109835(DDX5,p68 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109835(DDX5,p68 Antibody) at dilution of 1:50 (under 40x lens)