

# GPR81 antibody



Catalog Number: 111135

## Product name

GPR81 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR81 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human colon tissue.

## Preparation

This antibody was obtained by immunization of GPR81 recombinant protein (Accession Number: NM\_032554). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:111135(GPR81 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:111135(GPR81 Antibody) at dilution of 1:50 (under 40x lens)