

# neurturin-Specific antibody



Catalog Number: 113135

## Product name

neurturin-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

neurturin-Specific Rabbit Polyclonal antibody.  
Positive IHC detected in human colon tissue.  
Positive WB detected in human liver tissue, human heart tissue. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_004558).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113135(NRTN-Specific antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon using Catalog No:113135(NRTN-Specific antibody) at dilution of 1:50 (under 10x lens)