

PPM1G antibody



Catalog Number: 114135

Product name

PPM1G antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPM1G Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive IHC detected in human testis. Observed molecular weight by Western-blot: 70kd

Preparation

This antibody was obtained by immunization of PPM1G recombinant protein (Accession Number: NM_177983). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114135(PPM1G antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human testis using Catalog No:114135(PPM1G antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:114135(PPM1G antibody) at dilution of 1:100 (under 40x lens)