

TIMM8A antibody

Catalog Number: 116135

Product name

TIMM8A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TIMM8A Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human brain tissue, human liver tissue, mouse testis tissue, rat kidney tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 11 kDa

Preparation

This antibody was obtained by immunization of TIMM8A recombinant protein (Accession Number: NM_004085). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116135(TIMM8A antibody) at dilution of 1:500



IP Result of anti-TIMM8A (IP:Catalog No:116135, 4ug; Detection:Catalog No:116135 1:500) with mouse testis tissue lysate 4000ug.