

# HOXC10 antibody



Catalog Number: 111535

## Product name

HOXC10 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HOXC10 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, A549 cells. Observed molecular weight by Western-blot: 50kd

## Preparation

This antibody was obtained by immunization of HOXC10 recombinant protein (Accession Number: NM\_017409). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111535(HOXC10 antibody) at dilution of 1:500