

RAMP3 antibody

Catalog Number: 114535

Product name

RAMP3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RAMP3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse heart tissue, A549 cells, human heart tissue, human kidney tissue, human placenta tissue, Jurkat cells, MCF7 cells, mouse brain tissue, mouse lung tissue, mouse pancreas tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of RAMP3 recombinant protein (Accession Number: NM_005856). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114535(RAMP3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:114535(RAMP3 antibody) at dilution of 1:50 (under 10x lens)