

V5-tag antibody

Catalog Number: 107535

Product name

V5-tag antibody

Specificity

Human; other species not tested.

Antibody description

V5-tag Mouse Monoclonal antibody. Positive IP detected in Transfected HEK-293 cells. Positive WB detected in recombinant protein.

Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein G purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:5000-1:50000

IP: 1:2000-1:20000

Validations



Western blot of V5-tagged fusion protein with V5 tag monoclonal antibody (Catalog No:117343) at



various dilutions.

IP Result of anti-V5-tag (IP:Catalog No:117343, 4ug; Detection:Catalog No:117343 1:5000) with Transfected HEK-293 cells lysate 400ug.