

# Anti-GSK3 beta (3D10) Mouse antibody



Catalog Number: 168535

## Product name

Anti-GSK3 beta (3D10) Mouse antibody

## Specificity

Human, Mouse, Monkey, Rat

## Antibody description

GSK3 beta (3D10) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human GSK3B expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

## Validations



Western blot detection of GSK3 beta in 3T3, PC12 and A549 cell lysates using GSK3 beta mouse mAb (dilution 1:500). Predicted band size: 46KDa. Observed band size: 46KDa.