

# CPZ antibody



Catalog Number: 109535

## Product name

Polyclonal

CPZ antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

CPZ Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

## Dilutions

Recommended Dilution:

## Preparation


WB: 1:200-1:1000

This antibody was obtained by immunization of CPZ recombinant protein (Accession Number: BC006393). Purification method: Antigen affinity purified.

## Validations

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109535(CPZ Antibody) at dilution of 1:350

## Clonality