

# Anti-GRP78 BiP Rabbit antibody



Catalog Number: 168935

## Product name

Anti-GRP78 BiP Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GRP78 BiP Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human GRP78 BiP

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000-1/5000

IHC: 1/50-1/200

## Validations



Western blot detection of GRP78 BiP in C6 cell lysates using GRP78 BiP Rabbit mAb(1:1000 diluted).Predicted band size:72kDa.Observed band size:78kDa.