

Human UDP galactose-4'-epimerase / GALE (His Tag) recombinant protein



Catalog Number: 500235

General Information

Gene Name Synonym

Galactowaldenase; UDP-N-acetylgalactosamine 4-epimerase; UDP-N-acetylglucosamine 4-epimerase; UDP-galactose 4-epimerase

Protein Construction

A DNA sequence encoding the mature form of human GALE (Q14376) (Met1-Ala348) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human GALE consists of 363 amino acids and predicts a molecular mass of 40.1 KDa. It migrates as an approximately 36 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM HEPES, 150mM NaCl, 10% Glycerol, pH 7.5.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human GALE / UDP galactose-4-epimerase Protein (His Tag) SDS-PAGE