

EMB antibody

Catalog Number: 110235

Product name

EMB antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

EMB Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in HeLa cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of EMB recombinant protein (Accession Number: NM_198449). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110235(EMB Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110235(EMB Antibody) at dilution of 1:50 (under 40x lens)